


Dental Ratio
buy smart



coreflon

INTELLIGENT SUTURE



PTFE SURGICAL SUTURE

BLACK NEEDLE AND WHITE SUTURE				USP							
				3/0		4/0		5/0		6/0	
BLACK · DS Circle 1/2		Round	Length 55 cm	<input type="checkbox"/>		<input type="checkbox"/>	C5DS135512B	<input type="checkbox"/>	C6DS135512B*		
				<input type="checkbox"/>	C3DS165512B	<input type="checkbox"/>	C4DS165512B*	<input type="checkbox"/>	C5DS165512B	<input type="checkbox"/>	C6DS165512B*
				<input type="checkbox"/>	C3DS195512B	<input type="checkbox"/>	C4DS195512B*	<input type="checkbox"/>	C5DS195512B		
				<input type="checkbox"/>	C3DS225512B	<input type="checkbox"/>	C4DS225512B*				

SILVER NEEDLE AND WHITE SUTURE				USP							
				3/0		4/0		5/0		6/0	
SILVER · DS Circle 1/2		Round	Length 55 cm			<input type="checkbox"/>	C5DS135512S*	<input type="checkbox"/>	C6DS135512S*		
				<input type="checkbox"/>	C3DS165512S*	<input type="checkbox"/>	C4DS165512S	<input type="checkbox"/>	C5DS165512S*	<input type="checkbox"/>	C6DS165512S*
				<input type="checkbox"/>	C3DS195512S*	<input type="checkbox"/>	C4DS195512S	<input type="checkbox"/>	C5DS195512S*		
				<input type="checkbox"/>	C3DS225512S*	<input type="checkbox"/>	C4DS225512S				

SILVER · DM Circle 3/4		Flat	Length 55 cm			<input type="checkbox"/>	C5DM135512S	<input type="checkbox"/>	C6DM135512S		
				<input type="checkbox"/>	C3DM165512S*	<input type="checkbox"/>	C4DM165512S*	<input type="checkbox"/>	C5DM165512S	<input type="checkbox"/>	C6DM165512S*
				<input type="checkbox"/>	C3DM195512S*	<input type="checkbox"/>	C4DM195512S*	<input type="checkbox"/>	C5DM195512S*		
				<input type="checkbox"/>	C3DM225512S*	<input type="checkbox"/>	C4DM225512S*				

DS

Round needle with cutting blade - The needle is designed for surgical procedures where only minimal tissue traumatization is desired, but excellent penetration, better than with round needles, is still required. Provides a smooth transition between the hard tissues..

DM

Needle with reverse cutting blade type P - The needle has a specially designed tip, the third section of the length of the needle is sharp, while the remaining part has an oval cross-section, flattened for better grip in a vise.

* only on request available

DRS International GmbH · ISSUE 07/01/2014

	DRS International GmbH Raiffeisenstr. 15 a 40764 Langenfeld Germany Phone +49 2173 39800 0	info@dental-ratio.de www.dental-ratio.de Fax +49 2173 39800 20	Customer ID
	<input type="checkbox"/> Standard	<input type="checkbox"/> Express	Company
			Street
			Zip Code/City
			Phone
			E-Mail