

Quantity Art. No. Product Quantity Art. No. Product

**OKTAGON® TISSUE LEVEL · SOLID ABUTMENTS**

<input type="checkbox"/>	81038	TL RP Solid abutment yellow, 6°, H 4.0 mm	<input type="checkbox"/>	55283	TL RP Impression cap (snap-on)
<input type="checkbox"/>	81083	TL RP Solid abutment grey, 6°, H 5.5 mm	<input type="checkbox"/>	55288	TL RP Positioning cylinder yellow
<input type="checkbox"/>	81096	TL RP Solid abutment blue, 6°, H 7.0 mm	<input type="checkbox"/>	55299	TL RP Positioning cylinder grey
<input type="checkbox"/>	55478	TL RP Plastic coping for solid abutment, bridge	<input type="checkbox"/>	55394	TL RP Positioning cylinder blue
<input type="checkbox"/>	55487	TL RP Plastic coping for solid abutment, crown	<input type="checkbox"/>	55396	TL RP Solid abutment analog yellow
<input type="checkbox"/>	55668	TL RP Solid abutment driver 6°, L 15 mm, (S)	<input type="checkbox"/>	55397	TL RP Solid abutment analog grey
<input type="checkbox"/>	72088	TL RP Solid abutment driver 6°, L 21 mm, (M)	<input type="checkbox"/>	55400	TL RP Solid abutment analog blue

**OKTAGON® TISSUE LEVEL · BALL ATTACHMENTS**

<input type="checkbox"/>	38087	TL RP Abutment ball attachment GH 0.0 mm	<input type="checkbox"/>	82758	TL RP Matrix for ball attachment (incl. all O-Rings)
<input type="checkbox"/>	38071	TL RP Implant Analog for ball attachment	<input type="checkbox"/>	89611	TL RP O-Ring 70 shore red Pull-off weight <b>light</b>
<input type="checkbox"/>	38073	TL RP Ball attachment driver	<input type="checkbox"/>	89612	TL RP O-Ring 80 shore blue Pull-off weight <b>medium</b>
			<input type="checkbox"/>	89613	TL RP O-Ring 90 shore black Pull-off weight <b>strong</b>

Art. No. 90512 · DRS International GmbH · ISSUE October 2014



DRS International GmbH  
Raiffeisenstr. 15 a  
40764 Langenfeld  
Germany

info@dental-ratio.de  
www.dental-ratio.de

Phone +49 2173 39800 0

Fax +49 2173 39800 20

Standard

Express

Customer ID

Company

Name

Street

Zip Code/City

Phone

E-Mail